

MS Table 2B Hi-Flow

Номер части: 10-1028

Описание

Стол MS был разработан как эксклюзивный продукт Sciex для размещения систем Triple Quad™, QTRAP™ или IVD MS. Стол MS поставляется в стандартной и премиальной версиях, премиальная версия включает в себя встроенный Infinity 1031, позволяющий подавать завесу, исходный и выхлопной газы. Кроме того, MS Table также служит доступным звукопроницаемым корпусом для необходимых вакуумных насосов для этих приложений.




Generator Details

- **Типы газа:** Воздух
- **Встроенный компрессор:** Нет

Ключевая особенность

- Стол предназначен для размещения масс-спектрометра Sciex.
- Может поставляться отдельно или со встроенным бескомпрессорным генератором азота Infinity 1031.
- Уменьшение площади лаборатории
- Требуется минимальная настройка.
- Звукоизоляция, снижающая шум генератора.



Технические характеристики

MS Table 2B Hi-Flow

Тип газа	Воздух
Мин. давление на входе:	120psi/8,3bar
Максимальное давление на входе:	145psi/10bar
Частота:	50 / 60 Hz
Размер (ВхШхГ), мм	804 x 1000 x 830 mm
Размер (ВхШхГ), дюймы	31,7 x 39,4 x 32,7 inches

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

